


EXP. NUMBER 1	EXPERIMENT/SUBJECT Acid Base Titration	DATE 5/2/11	01
NAME Dr. Grant EDWARDS	LAB PARTNER None	LOCKER/DESK NO. 22	COURSE & SECTION NO. 120 L - 16

Unknown # 25

Objective: To titrate an unknown sample of acid with a standardized base solution and determine the CONCENTRATION of the unknown in Molarity.



### PROCEDURE:

- Obtain a solution of standardized base.
- Prime a buret, ~~by~~ with ~10 mL of Base
- Acquire a sample and record unknown number (also at top!)
- Quantitatively transfer unknown to 100 mL volumetric flask and fill to the line with distilled water.
- ~~USE A pipet to Transfer~~ prime a ~~10.00~~ <sup>15.00 mL</sup> mL pipet with the unknown Acid solution
- pipet <sup>15.00 mL</sup> 10.00 mL of Acid into an Erlenmeyer flask, Add ~2 drops phenolphthalein indicator, and ~25 mL distilled water.
- Titrate the erlenmeyer flask mixture with the base in the buret until the solution is very light pink
- Repeat until reproducible (for more trials).

### Observations:

Concentration  $\text{NaOH}$  0.1023 M

unknown number 25

Better: also record observations. Solution colors, smell, etc? Especially if changes occur!

Changed per prelab instruction

Note changes to procedure as needed.

Trial data for each trial in table on next page.

SIGNATURE	DATE	WITNESS/TA	DATE
-----------	------	------------	------

EXP. NUMBER	EXPERIMENT/SUBJECT	DATE	02	
NAME	LAB PARTNER	LOCKER/DESK NO.	COURSE & SECTION NO.	

Data:

Table 1: Titration Trials

Record data, be clear.

Trial #	T1	T2	T3	T4	T5
Amount of Acid (mL)	<del>10.00</del> 15.00	15.00	15.00	15.00	15.00
Initial buret reading (mL)	0.00	13.24	0.05	12.06	24.06
Final buret reading (mL)	13.24	25.74	12.06	24.06	36.09
Volume base (mL)	13.24	12.50	12.01	12.00	12.03 mL
Molarity of Unknown Acid (M)			0.08225M	0.08184	0.08204
Endpoint Color	dark pink	dark pink	light pink	light pink	light pink
Valid or Not, why.	No past endpoint	No	good	good	good

buret dripped to fast!

Shown calculations are clear with units and easy to follow.

Calculations

T3 Molarity of unknown Acid

$$\frac{0.1023 \text{ mol NaOH}}{1 \text{ L sol'n}} \times \frac{12.06 \text{ mL}}{13.24 \text{ mL sol'n}} \times \frac{1000 \text{ mL sol'n}}{1 \text{ L sol'n}}$$

$$\frac{1 \text{ mol Acetic Acid}}{1 \text{ mol NaOH}} \times \frac{15.00 \text{ mL AA sol'n}}{1000 \text{ mL}} = 0.0822492 \frac{\text{mol AA}}{\text{L sol'n}} \Rightarrow \boxed{0.08225 \text{ M AA}}$$

You can also have calculation result in tables. Replicate trials may be calculated using excel, but show calculation in notebook also to check excel formula results.

Trials 3 → 4 standard deviations

- T1 0.08225
- T2 0.08184
- T3 0.08204

$$\bar{x} = 0.082043 \Rightarrow \boxed{0.08204 = \text{Average}}$$

Std deviation

$$s = \sqrt{\frac{\sum(x-\bar{x})^2}{n-1}} = 0.0002050$$

$$\% \text{RSD} = \frac{0.0002050}{0.08204} \times 100 = \boxed{0.2487\%}$$

These calculations completed in Microsoft Excel.

You don't have to show average or standard deviation calculations as I assume you can do them correctly. However you must be clear about what trials are being used and why.

SIGNATURE	DATE	WITNESS/TA	DATE
-----------	------	------------	------

EXP. NUMBER	EXPERIMENT/SUBJECT	DATE	03
NAME	LAB PARTNER	LOCKER/DESK NO.	COURSE & SECTION NO.

Calculations Continued.

95% Confidence Interval.

$$\bar{x} \pm \frac{t \cdot s}{\sqrt{n}} \Rightarrow 0.08204 \text{ M} \pm \frac{4.303 \cdot 0.0002050}{\sqrt{3}} \Rightarrow 0.08204 \text{ M} \pm 0.00051 \text{ M}$$

t from table of Students t values for 95% confidence from text book.

Results.

Table 1 contains the molarity of the unknown Acid for Trials 3, 4, 5. The average molarity was found to be 0.08204 M with a standard deviation of 0.0002050 M. The 95% confidence interval for these three trials was 0.08204 M  $\pm$  0.00051 M.

Conclusions.

For Sample # 25 the Average molarity was found to be 0.08204 M. Acetic Acid with a standard deviation of 0.0002050 M over 3 trials which has a 95% confidence limit of 0.08204 M  $\pm$  0.00051 M.

Discussion:

In this experiment, a standardized solution of NaOH was used to titrate an unknown solution of acetic acid. The 95% confidence interval of 0.08204 M  $\pm$  0.00051 M says that with 95% confidence, the real concentration of the acid is between 0.08153 M and 0.08255 M.

One possible error that could have occurred was with the given concentration of the NaOH, if this was incorrect the resulting molarity found would also be correct. a better method would be to experimentally standardize the solution with a primary standard acid.


The second source of error could be the 15.00 mL pipet and/or the buret. Volume readings measurements are very important for both the unknown acid and standardized base.

SIGNATURE	DATE	WITNESS/TA	DATE
-----------	------	------------	------

EXP. NUMBER	EXPERIMENT/SUBJECT	DATE	04	
NAME	LAB PARTNER	LOCKER/DESK NO.	COURSE & SECTION NO.	

Discussion Contd.

A third source of error was the indicator used. If the indicator changed colors after the endpoint, actual equivalence point, the calculated molarity of the Acid would be too high. This is why the proper selection of an Indicator is key.



SIGNATURE	DATE	WITNESS/TA	DATE
-----------	------	------------	------